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Technical Notes

Collagenase NB 6 GMP Grade For Isolation and Passaging of Stem Cells

Collagenase NB 6 GMP Grade for Isolation and Passaging of Stem Cells

Collagenase plays a crucial role in isolation and passaging of stem cells dedicated for transplantation into humans. Nordmark Arzneimittel provides **Collagenase NB 6 GMP Grade** particularly for these clinical applications. This enzyme integrates superior quality with easy handling. It is suitable for isolation of a broad variety of cells, including stem cells (e.g. ADSC), and for stem cell passaging (e.g. hESC).

Collagenase NB 6 GMP Grade is sterile thus it is ready for use. It contains collagenase classes I and II as well as proteolytic activities such as neutral protease and clostripain. Therefore, Collagenase NB 6 GMP Grade is a mild and effective enzyme producing high yields of viable cells.

GMP Compliant Manufacturing

Isolation and passaging of cells dedicated for tissue engineering and transplantation into humans require a collagenase with reliable quality. For this reason Collagenase NB 6 GMP Grade is manufactured in compliance with the EU guide

to Good Manufacturing Practice (GMP) for active pharmaceutical ingredient by the German pharmaceutical company Nordmark.

Safety

Collagenase NB 6 GMP Grade meets high safety standards since TSE safety of the manufacturing process is certified by the EDQM. Each lot possesses low endotoxin level and is tested for abnormal toxicity according to European Pharmacopoeia. In addition, data for virus validation and stability studies according to ICH guidelines are available.

Collagenase NB 6 GMP Grade – The superior collagenase for clinical applications

Technical Notes

Collagenase AF-1 GMP Grade Neutral Protease AF GMP Grade Safe, Gentle & Free of Animal-Derived Materials

Animal-free GMP Grade Collagenase and Neutral Protease for Cell Isolation

Nordmark, in close collaboration with the company SERVA, has developed the world's first highly purified collagenase and purified neutral protease **completely free of any animal-based components**. Using raw materials from animal origin often causes concerns regarding the introduction of potential animal-derived pathogens. To rule out these concerns, **Collagenase AF-1 GMP Grade** and **Neutral Protease AF GMP Grade** are derived from *Clostridium histolyticum* grown in a medium containing carefully selected **plant-based** ingredients. Furthermore, the sophisticated production process ensures that there is no cross-contamination with any animal-derived materials.

Therefore, these enzymes provide the **highest level of safety and reduced regulatory hurdles**.

Highly Purified for Exact Dosage and High Enzyme Stability

Collagenase AF-1 GMP Grade and Neutral Protease AF GMP Grade are highly active and defined products designed to isolate **especially sensitive cells** from tissues such as liver, dental pulp or pancreas.

Collagenase AF-1 GMP Grade is chromatographically highly purified containing:

- > Class I and Class II collagenase
- > Very low levels of endotoxins
- > Reduced proteolytic side activities

Often a certain amount of additional enzymatic side activities is essential. For this reason, chromatographically purified Neutral Protease AF GMP Grade is provided separately, which has the great benefit of:

- > Individual and precise dosage
- > Higher stability of the enzymes
- > Optimized yield and viability of sensitive cells

Tested for Islet Isolation by Leading Scientists

The efficacy and reliable performance of these **new animal-free enzymes** have been **confirmed by leading experts** in the field of islet isolation from renowned institutions such as University of Geneva, University of California, San Francisco, and Uppsala University.

Islets isolated with AF enzymes. Photo: Engbert, EGP

Manufactured According to GMP Guidelines

Our animal-free enzymes are produced with the expertise of the world's largest manufacturer of collagenase, the pharmaceutical company Nordmark. Compliance with GMP guidelines guarantees the **highest quality** and a reliable **lot-to-lot consistency**.

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© 2004 Blackwell Publishing Ltd, *Journal of Internal Medicine* 255: 111–118